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


Automate data analysis and simplify MS<sup>n</sup> data analysis through customizable workflows. An example of a customizable workflow with eight separate searches on the same raw file.

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

## • Volcano plots (New in PD 2.2)

## • Pathway maps (New in PD 2.2)

Proteome Discoverer Tools for statistical and biological interpretation. Includes study management for mapping of quantitation channels to biological information about the sample, with fast, robust algorithms for peak selection and peptide quantification to extract the maximum possible information from quantitative proteomics studies.

Third party nodes for data processing, facilitates the processing of large datasets. Additionally, biological context by automatically annotating identified proteins with GO classifications, PTMs sites, and literature references from public databases (active maintenance required) can be illustrated.

Effective location of low-level sample components using creative data interrogation strategies such as generating inclusion and exclusion lists and importing those lists directly into instrument methods is provided, and data comparison, exchange, and verification with support for data standards developed by the HUPO Proteomics Standards Initiative (PSI) is additionally supported.

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**Taxonomy** [ ] Primates

**Enzyme** TrypsinP      **Allow up to** 2 missed cleavages

**Quantitation** None

**Crosslinking** Disulfide bridge in Elongation factor 2 in human

**Fixed modifications** Carbamidomethyl (C)

Display all modifications

**Variable modifications** Oxidation (M)

**Peptide tol.** 1 ppm @ 1<sup>st</sup> C

**MS/MS tol.** 1 ppm

**Peptide charge** 2+      **Monoisotopic** ☒ **Average** ☐

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**Data format** Mascot generic      **Precursor** ☐ m/z

**Instrument** CID+ETD      **Error tolerant** ☐

**Decoy** ☐ **Report top** AUTO hits

**Start Search** **Reset Form**

[illegible]

[illegible]

- 
- The screenshot shows a Windows desktop with a taskbar at the bottom containing icons for Internet Explorer, Google Chrome, and several folders. The active window is a web browser displaying a search engine results page for the query "IP Address". The search results include a link to "IP Address - Wikipedia" and a table of IP addresses.
- The table has the following columns: **IP Address**, **OS**, **Resolving**, **Physical Memory**, **Swap Size**, and **Disk Space**. The table contains 10 rows of data, each representing a different IP address and its associated system information.
- | IP Address   | OS    | Resolving | Physical Memory | Swap Size | Disk Space |
|--------------|-------|-----------|-----------------|-----------|------------|
| 192.168.1.1  | Linux | 0.0       | 4.75 GB         | 0.00 GB   | 4.80 GB    |
| 192.168.1.2  | Linux | 0         | 4.75 GB         | 0.00 GB   | 4.80 GB    |
| 192.168.1.3  | Linux | 0         | 4.75 GB         | 0.00 GB   | 4.80 GB    |
| 192.168.1.4  | Linux | 0         | 4.75 GB         | 0.00 GB   | 4.80 GB    |
| 192.168.1.5  | Linux | 0         | 4.75 GB         | 0.00 GB   | 4.80 GB    |
| 192.168.1.6  | Linux | 0         | 4.75 GB         | 0.00 GB   | 4.80 GB    |
| 192.168.1.7  | Linux | 0         | 4.75 GB         | 0.00 GB   | 4.80 GB    |
| 192.168.1.8  | Linux | 0         | 4.75 GB         | 0.00 GB   | 4.80 GB    |
| 192.168.1.9  | Linux | 0         | 4.75 GB         | 0.00 GB   | 4.80 GB    |
| 192.168.1.10 | Linux | 0         | 4.75 GB         | 0.00 GB   | 4.80 GB    |

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| #  | Visible?                            | Starred?                            | Bio View:<br>Identified Proteins (32) | Accession N | Molecular We | Protein Group | Taxonomy | biological reg | cellular proce | development? | metabolic pro | multicellular | response to s | cytoplasm | cytoskeleton | endosome | extracellular | membrane | mitochondriol | nucleus | organelle par | plasma mem | binding | catalytic acti | enzyme regul | molecular fun | structural mc | 1_M | 1_G |
|----|-------------------------------------|-------------------------------------|---------------------------------------|-------------|--------------|---------------|----------|----------------|----------------|--------------|---------------|---------------|---------------|-----------|--------------|----------|---------------|----------|---------------|---------|---------------|------------|---------|----------------|--------------|---------------|---------------|-----|-----|
| 1  | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | (P02470) Alpha crystallin A ch...     | CRYAA_BOVIN | 20 kDa       | Bos taurus    |          |                |                |              |               |               |               |           |              |          |               |          |               |         |               |            |         |                |              |               | 27            | 29  |     |
| 2  | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | (P07318) Beta crystallin B1           | CRBB1_BOVIN | 28 kDa       | Bos taurus    |          |                |                |              |               |               |               |           |              |          |               |          |               |         |               |            |         |                |              |               | 32            | 33  |     |
| 3  | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | (P02470) Alpha crystallin A ch...     | CRYAA_BOVIN | 20 kDa       | Bos taurus    |          |                |                |              |               |               |               |           |              |          |               |          |               |         |               |            |         |                |              |               | 27            | 29  |     |
| 4  | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | (P07318) Beta crystallin B1           | CRBB1_BOVIN | 28 kDa       | Bos taurus    |          |                |                |              |               |               |               |           |              |          |               |          |               |         |               |            |         |                |              |               | 32            | 33  |     |
| 5  | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | (P02470) Alpha crystallin A ch...     | CRYAA_BOVIN | 20 kDa       | Bos taurus    |          |                |                |              |               |               |               |           |              |          |               |          |               |         |               |            |         |                |              |               | 27            | 29  |     |
| 6  | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | (P07318) Beta crystallin B1           | CRBB1_BOVIN | 28 kDa       | Bos taurus    |          |                |                |              |               |               |               |           |              |          |               |          |               |         |               |            |         |                |              |               | 32            | 33  |     |
| 7  | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | (P02470) Alpha crystallin A ch...     | CRYAA_BOVIN | 20 kDa       | Bos taurus    |          |                |                |              |               |               |               |           |              |          |               |          |               |         |               |            |         |                |              |               | 27            | 29  |     |
| 8  | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | (P07318) Beta crystallin B1           | CRBB1_BOVIN | 28 kDa       | Bos taurus    |          |                |                |              |               |               |               |           |              |          |               |          |               |         |               |            |         |                |              |               | 32            | 33  |     |
| 9  | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | (P02470) Alpha crystallin A ch...     | CRYAA_BOVIN | 20 kDa       | Bos taurus    |          |                |                |              |               |               |               |           |              |          |               |          |               |         |               |            |         |                |              |               | 27            | 29  |     |
| 10 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | (P07318) Beta crystallin B1           | CRBB1_BOVIN | 28 kDa       | Bos taurus    |          |                |                |              |               |               |               |           |              |          |               |          |               |         |               |            |         |                |              |               | 32            | 33  |     |

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The screenshot shows the KEGG Pathway database interface. The search bar contains 'Fabp1' and the results show '15 pathways found'. A red arrow points from the text 'Pathway 情報サイトへのリンク' to the search bar.

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